

Procedure for Rietveld refinement

Atsushi Kubo

Selection of proper 2D data

- Check Debye rings. It should be less spotty and continuous.
- Check cake pattern. It should not be wavy.

Data conversion of 2D data to 1D data using Fit2D

- Mask spots that have saturated or extremely strong intensity.
- Mask 2D area to fan-shape for integration to maximize 2θ range in 1D data.

Correction of diffraction intensity by taking x-ray absorption by diamond into account (using abs_coef_carbon.xls)

Determine I_0 at each data point based on the following equation.

$$I_{\text{obs}} = I_0 e^{-\mu x}$$

- μ ; a function of X - ray energy and density of the material (in this case diamond)
- x ; thickness of the material

Determination of Estimated Standard Deviation (ESD) of intensity at each 2θ data point (using abs_coef_carbon.xls)

Obtain σ at each data point using the equation suggested by Von Dreele (CCP14):

$$\sigma = \frac{F\sqrt{I}}{\tan 2\theta}$$

σ : ESD

F: constant

Convert 3-columns text data in to text file for GSAS ESD format (using Powder4)

Load data and add phase(s) into GSAS/EXPGUI

- Add phase(s) with space group and initial lattice parameters
- Add atom(s) and input initial fractional coordinates of the atom(s)
- Load the data file of (2θ , I_0 , σ) (GSAS ESD format)
- Load instrument parameter file (e.g. polarization factor of x-ray = 0.95)

Start Le Bail refinement

- Set data limits and excluded regions
- Fit background graphically (using typically 10 terms. Background is typically difficult to refine for angle-dispersive data!!)
- Refine profile terms (= curve shape parameters in Pseudo-Voigt function)
 - Refine only GW and LX. Start from GW, then refine both GW and LX.

$$\sigma^2 = U \tan^2 \theta + V \tan \theta + W$$

σ^2 : Variance of the peak in pseudo - Voigt function

(Here $U = GU$, $V = GV$, $W = GW$)

$$\gamma = \frac{X}{\cos \theta} + Y \tan \theta$$

γ : Lorentzian coefficient in pseudo - Voigt function

(Here $X = LX$, $Y = LY$)

Refine lattice parameters

Start Rietveld refinement

(1-i) fix lattice parameters and profile terms, refine phase fractions.

(1-ii) fix lattice parameters and profile terms, refine both phase fractions and atomic positions.

(1-iii) fix lattice parameters and profile terms, refine phase fractions, atomic positions, and displacement parameters U_{iso} simultaneously (but refinement of displacement parameters tend to yield physically meaningless values, so typically we need to fix U_{iso} to reasonable values and never refine).

(2) fix all the above parameters, and refine spherical harmonic terms for preferred orientation correction.

(3) refine all the above parameters together.